## 3. Multi-channel pipette - all models

Type: Basic (body)

Quick notes: · Shaft is supplied as a kit, including compression spring & washer (that guarantee back compatibility) • Since 2014/2015 no external shafts are used, replaced with new version of shafts, springs and shaft box (black)

Item	Name
Α	Pushbutton (cap & knob)
A3	Pipeting arbor w/bush
for parts (B) check the 1ch pipettes	
7	Manifold nut (set)
8	Piston
9	Piston spring
16	O-ring
17	Shaft kit***
19	Tips ejector
19A	Ejector pusher
23	Manifold housing set
Ν	Ejector button
N1	Ejector arbor
N2	Ejector arbor fixing kit
Optional:	
18	Shafts assembly
24	Manifold housing w/ejector
34H	Screws for housing set*
34E	Ejector mounting set**
34P	Screw M2.2×9.5
34S	Screw M2.2×13
Κ	Calibration key
Item code matches the character used	

in service manual, spare parts list and user manual

\* set = 5pcs of each screw size: 4/6/10/16mm

\*\* set = 1pc of each element:

- top screw+washer

- bottom screw+spring+holder

