


Experience the Picus[®] 2 Difference


How Picus[®] 2 Benefits Your Pipetting



Connectivity

Use Sartorius pipetting mobile app to smoothly run sample preparation workflows and adjust the pipette setting automatically, taking your productivity to the next level.

Open connectivity enables you to integrate Picus[®] 2 to be part of an existing system.




2X Faster Than Mechanical Pipettes

Thanks to multi-dispensing mode and other modes work steps are reduced for both single and multichannel models



8 Modes for Versatile and Demanding Pipetting Applications

Revolutionize repetitive work and utilize multiple pipetting modes, from pipetting and diluting to serial dispensing to speed up work—compared to only two modes for the mechanical pipettes.



Best Pipetting Results Independent of User Experience:

Fully electronic pipetting functions and more precise and smooth piston movement compared to mechanical pipettes.

How Picus[®] 2 Benefits You



100g

Weight makes Picus[®] 2 as Light as Mechanical Pipettes

Light weight and ergonomic handle of Picus[®] 2 makes it easy to use and helps reduce fatigue during pipetting, especially in long pipetting series.



95%

Lighter Operational Forces Help Reduce Hand and Wrist Strain and Injuries

Fully electronic operation of Picus[®] 2 eliminates most of the strain that occurs due to the force applied by the thumb for pipetting and tip ejection in mechanical pipettes.



1:1,000

One Hour Battery Charge 1,000 Pipetting Cycles

Always ready for pipetting as fast charging enables convenient pipetting throughout the day.