


**Spare parts
Ersatzteile
Pieces detachees**

1062800 1µl-10ml



1



4

1062930 MCP





2*



3*

1-10 ml

0.5-5 ml




- * 2. Cap
- * 3. Handle
- * 12. O-ring
- * 42. Tip Cone assy


All 0.2µl-10ml

- 2. 2214550
- 3. 2214650 **1-10ml**
- 3. 2214640 **0,5-5ml**
- 12. 1033050 **10ml**
- 12. 1030230 **5ml**
- 42. 2214540 **1-10ml**
- 42. 2214530 **0.5-5ml**


100-1000µl



20-200µl



10-100µl



- * 2. Cap
- * 3. Handle
- * 17. O-ring
- * 42. Tip Cone assy

100-1000µl

- 3. 2214630
- 17. 1030020
- 42. 2214520


20-200µl

- 3. 2214620
- 17. 1030160
- 42. 2214510


10-100µl

- 3. 2214610
- 17. 1030510
- 42. 2214500

2-20µl / 5-50µl



1-10µl



- * 2. Cap
- * 3. Handle
- * 19. O-ring
- * 20. O-ring
- * 24. O-ring
- * 42. Tip Cone assy

5-50µl

- 3. 2214600
- 19. 1030500
- 20. 1033060
- 42. 2214490


2-20µl

- 3. 2214660
- 19. 1030380
- 20. 1033110
- 42. 2214480


1-10µl

- 3. 2214560
- 19. 1030380
- 20. 1030060
- 24. 1030170
- 42. 2214470

0.5-5 µl



0.2-2 µl



- * 2. Cap
- * 3. Handle
- * 19. O-ring
- * 20. O-ring
- * 24. O-ring
- * 28. Sealing
Combination
- * 42. Tip Cone assy




0.5-5 µl

- 3. 2214590
- 19. 1030380
- 20. 1033160
- 24. 1030170
- 42. 2214460

0.2-2 µl

- 3. 2214580
- 24. 1030170
- 28. 2205730
- 42. 2214450

*** 50**



31

33

32

35

36

37

38

39

40

41

30

42*

1-10µl


- 2. 2214550
- 3. 2214560
- 42. 2209110 **12 pcs**
- 50. 2214670 **8-ch**
- 50. 2214710 **12-ch**



23

24

*** 50**



*** 2. Cap**
*** 3. Handle**
*** 42. Tip Cone assy**
*** 50. Module**





16



17

18

31



20

21



33

32

34

35

36

37

30



42*

23



5-50µl

2. 2214550

3. 2214600

42. 2212750 **12 pcs**

50. 2214680 **8-ch**

50. 2214720 **12-ch**

*** 50**



*** 2. Cap**


*** 3. Handle**

*** 42. Tip Cone assy**

*** 50. Module**



42*



10-100µl

2. 2214550


3. 2214610

42. 2212440 **12 pcs**

50. 2214690 **8-ch**

50. 2214730 **12-ch**

*** 50**



13

16

17

31

33

32

34

35

36

30

42*

*** 2. Cap**

*** 3. Handle**

*** 42. Tip Cone Assy**

*** 50. Module**



18

20

21

22

23

30-300µl

2. 2214550

3. 2214750

42. 2209350 **12 pcs**

50. 2214700 **8-ch**

50. 2214740 **12-ch**